



- 1 Entry port Ø 50 mm \*
- 2 Not assigned
- 3 Fresh air
- 4 C-Terminal
- 5 Interface RS 232

- 6 Potential free contact
- 7 Main switch
- 8 Adjustable operating temperature limiter\*
- 9 Pushbutton for lamp frame\*
- 10 Switch for test space illumination
- 11 Connector panel\*
- 12 Connection for automatic water replenishment
- 13 Drain for test space and humidification water
- 15 Mains supply, cable length approx. 5 m
- 16 Door with window
- 17 Irradiation module 2\* and 3\*
- 18 Connection for irradiation module 2\* and 3\*

For narrow entry openings components (e.g. hinges) may be screwed off to minimize width.